

CoVent Prime



ISO 9001 : 2015
ISO 13485:2016

Panakeia Medtech Pvt. Ltd.

AMTZ Campus, VM Steel Projects, Pragathi Maidan,

Visakhapatnam 530 031 (AP)

Mob: 9811340469 E-mail: panakeia.india@gmail.com

CoVent Prime

FEATURE

- 12.1" TFT color screen
 - Easy setting by navigator wheel knob and touch key
 - Wide choice of ventilation modes including CPAP, BILEVEL
 - Flow and pressure trigger
 - Static compliance and resistance monitor
 - PEEP, Pressure - Time, Volume - Time, Flow - Time
 - Three - level visual and audible alarm, with concise words describe
 - Advanced built in electronic air -oxygen mixing device
 - Durable and accurate built-in flow sensor
 - integrated expiration valve,
 - Backup apnea ventilation
 - Inbuilt nebulizer
 - Built in battery
- CE & ISO certified



CoVent Prime

CMV / VCV MODE



In VCV mode, the ventilator delivers mandatory breaths by setting tidal volume. When the ventilator detects patient's inspiratory effort, it delivers a patient-initiated mandatory (PIM) breathe. If the ventilator does not detect inspiratory effort, it delivers a ventilator-initiated mandatory (VIM) breath at an interval based on the set respiratory rate. Breaths can be pressure-triggered or flow - triggered in VCV mode.

PSV MODE

In PSV mode, inspiration is usually initiated by patient effort. Breaths are initiated via pressure of flow triggering. Clinician can also initiate a manual inspiration during PSV. PSV breaths are Controlled by pressure (preset PSV level + PEEP)
Limited by pressure (preset PSV LEVEL + PEEP + margin)
Cycled by time (PSV Tmax) or flow (PSV Cycle)



SIMV MODE



SIMV is a mixed ventilation mode that allows both mandatory and spontaneous breaths. The Mandatory breaths can be volume - based (SIMV-VC) or pressure-based (SIMV-PC), and the spontaneous breaths can be pressure-assisted. clinician can select pressure-triggering or Volume-triggering in SIMV

CoVent Prime

The 'Panakeia' 'Covent Prime' is electronically controlled and pneumatically operated ventilator for pediatric and adult patients.

It has different volume and pressure control modes apart from support modes like Pressure and Volume support.

It also provides for Pressure Regulated Volume Adjust.

All the set and patient parameters are displayed in LCD screen.

The monitored data include expired tidal volume, expired minute volume, PIP, MAP, Plateau Pressure and delivered FiO₂. The display also provides for pressure, volume and flow curve in real time.

The ventilator also provides for Nebulization and comprehensive alarm systems.

A low noise air compressor has been integrated in the unit.

SPECIFICATIONS

Gas supply	: Oxygen and air at 3 to 4 bar
Modes	: CMV - VC, CMV -PC, ACV, PS, 2PAP+PS, Standby SIMV (VC), SIMV (PC) + PS,
FiO ₂	: 21 to 100%
Tidal volume	: 20 to 1500 ml
Frequency	: 1 to 60 BPM
I:E Ratio	: 1:4 to 4:1
Inspiratory time	: 0.3 to 3 seconds
Inspiratory hold	: 0 to 2 seconds
Pick flow	: 150 Liters/min
PEEP	: 0 to 30 cmH ₂ O
Sensitivity	: 1 -10 cmH ₂ O (Pressure), 1-10 LPM (Flow)
PSV	: 0 – 30 cmH ₂ O
2PAP EPAP	: 0-30,
IPAP	: 0-40cmH ₂ O
Alarms	: Oxygen/Air Failure, Patient Disconnection
High Pressure	
Monitored Patient Data	: Exp. Tidal Volume, PIP, Respiratory Rate MAP, FiO ₂
Curves Pressure	: Pressure Time, Flow-Time & Volume-Time
Loops	: Pressure - Volume, Flow - Volume
Electrical requirement	: 200 – 220 VAC, 50 Hz
Battery Back Up	: Two Hours on full charged condition
Accessories Standard	: Breathing circuit, Heated Humidifier, Oxygen Sensor, Air Compressor. Pressure Hose and Moisture Trap.

