

ANAESTHESIA MACHINE

ANAESTHESIA MACHINE (I)



Anti-Hypoxic pneumatic device -
Ensures no nitrous oxide flow if oxygen
Flow is lower than 0.5 ltr/min approx.
Minimum 25% (+ - %5 oxygen) in gas
mixture gas delivery.

Ratio control :1 :3 mixed ratio between
O₂. & N₂O gas Flow

Rotameters:" ISI & ISO Standard

Oxygen driven N₂O supply.

Audio visual Oxygen failure protection
device {OFPD}.

Higher efficiency forged regulators &
yokes with stainless steel fittings High
precision pressures relief and non-
return valve.

Goldman Halothane Vaporizers" Change
over mechanism 'with oxygen flush
(open & closed circuit).

Gas inlet Port O₂ & N₂O each 2 nos. Pin-
indexed gas cylinder Connections & N₂O
each 2nos & diameter index gas specific
Central Pipeline Connection each one for
O₂. N₂O



ANAESTHESIA MACHINE (II)

MANUFACTURED BY

Panakeia Medtech Pvt. Ltd.



Standard Bain's circuit for adult.

Trolley made of M.S. multi colour Powder coated finished, Stainless steel tabletop tray.

Modified monitoring accessories tray, drawer Top stainless steel, tray.

5" Dia Antistatic castor wheel with break

Oxygen outlet for driving ventilator.



ANAESTHESIA MACHINE (III)

Large Dia (63 mm) Panel Mounted cylinder pressure, operating pressure & pipeline pressure display Gauges for O₂ & N₂O.

Upgrade ability for: Flow Tec Vaporizer, Circle Absorber & Ventilator {Optional at extra cost}.

Size 1585mmX850mmX735mm (H x W x D)

Address:

Panakeia Medtech Pvt. Ltd.

AMTZ Campus, VM Steel Projects,

Pragathi Maidan,

Visakhapatnam 530 031 (AP)

Ph: 9811340469